


## Faeces Sample for Culture

### THIS IS WHAT YOU NEED


Testing kit given to you by your care provider. Collection paper or a thoroughly cleaned vessel.

### HOW TO DO IT



Refrain from wiping yourself after having a bowel movement, or collect faeces on a paper/in a thoroughly cleaned vessel. Adhere to the sample-taking intervals prescribed by your care provider.



Take out the tube and the swab from the packaging. Unscrew the tube's lid. The tube contains a non-dangerous liquid. The liquid is to remain in the tube.




Swipe the swab in the anus right after having a bowel movement, or dip the swab in collected faeces. There should be visible faeces on the swab.




Put the swab into the tube.

Break the swab at the indentation.

Carefully screw the lid back on.



Label the tube with your full personal identity number and name, and the date and time at which the sample was taken. Use the label given to you by your care provider.



Hand in the sample to your care provider on the same day as the sample was taken at the latest. Transport the sample in a plastic bag of ice. Bring valid ID.