

Missbruksanalyser i urinprov

Gränsvärden (cut-off nivåer)

SCREENING (screen) (immunokemisk analys)	Gränsvärde	VERIFIKATION (verif) (masspektrometrisk analys)	Gränsvärde
U-Amfetaminer/Ecstasy	500 µg/L	U-Amfetaminer Amfetamin Metamfetamin MDMA MDA	200 µg/L 200 µg/L 200 µg/L 200 µg/L
U-Bensodiazepiner	200 µg/L	U-Bensodiazepiner Desmetyldiazepam Temazepam Oxazepam Lorazepam 7-Amino-klonazepam 7-Amino-nitrazepam 7-Amino-flunitrazepam alfa-OH-midazolam alfa-OH-alprazolam	50 µg/L 50 µg/L 50 µg/L 50 µg/L 50 µg/L 50 µg/L 50 µg/L 50 µg/L 50 µg/L
U-Buprenorfin	5 µg/L	U-Buprenorfin Buprenorfin Norbuprenorfin	5 µg/L 5 µg/L
U-Cannabis	20 µg/L	U-Cannabis THC-COOH	10 µg/L
U-Dextropropoxyfen	300 µg/L	U-Dextropropoxyfen Nordextropropoxyfen	150 µg/L
U-Etylglukuronid (EtG)	0.5 mg/L	U-Etylglukuronid (EtG) Etylglukuronid (EtG) Etylsulfat (EtS)	0.5 mg/L 0.1 mg/L
U-Fentanyl	1 µg/L	U-Fentanyl Fentanyl Norfentanyl	1 µg/L 2 µg/L
U-GHB	10 mg/L	U-GHB	10 mg/L
U-Kokain	150 µg/L	U-Kokain Bensoylekgonin	100 µg/L
U-LSD	1 µg/L	U-LSD LSD 2-oxo-3-OH-LSD	0.3 µg/L 0.1 µg/L
U-Metadon	300 µg/L	U-Metadon Metadon EDDP	150 µg/L 150 µg/L
U-Opiater (total)	300 µg/L	U-Opiater (total) Morfin Kodein Etylmorfin	300 µg/L 300 µg/L 50 µg/L

U-Oxikodon	300 µg/L	U-Oxikodon	
		Oxikodon	150 µg/L
U-Tramadol	100 µg/L	U-Tramadol	
		Tramadol	50 µg/L
		O-desmetyltramadol	50 µg/L
<i>Endast masspektrometrisk paketanalys</i>			
		U-Missbruksläkemedel	
		Zolpidem	50 µg/L
		Zolpidem-4-syra	50 µg/L
		Zopiklon	50 µg/L
		Zaleplon	50 µg/L
		Pregabalin	100 µg/L
		Gabapentin	100 µg/L
		Metylfenidat	100 µg/L
		Ritalinsyra	100 µg/L
		Ketamin	50 µg/L
		Tapentadol	100 µg/L